

Jazzy[®] SELECT

The highly manoeuvrable Jazzy[®] Select delivers the advantages of in-line motor technology for enhanced efficiency, torque, range, and performance. Plus, it's loaded with performance and convenience features, like standard Active-Trac[®] Suspension, that make it as easy to enjoy as it is to use.



Jazzy[®] Select
shown in Red

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL
(Tel) 01869 324 600 • (Fax) 01869 323 070 • www.pride-mobility.co.uk

BUILT WITH PRIDE
Pride[®]
Mobility Products Ltd.

SPECIFICATIONS

Weight capacity	136 kg (21 Stone)
Maximum speed	Up to 6 km/h (4 mph)
Ground clearance	4.75 cm (1.875") at midframe
Turning radius	51 cm (20")
Overall length	90 cm (35.5") without front riggings
Overall width	60 cm (23.75")
Wheels	Drive: 25 cm (10")
Front anti-tips	13 cm (5") solid
Rear castors	15 cm (6") solid
Suspension	Active-Trac®
Drive train	Two-motor, in-line, mid-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 2 controller
Specialty controls	N/A
Battery charger	Off-board, 2 amp
Per-charge range	Up to 24 km (Up to 15 miles)
Battery requirements	(2) 12 volt deep cycle
Battery size	30 amp
Battery weight	11 kg (24.5 lbs.) each
Base weight	45 kg (99 lbs.)
Standard seat weight	17 kg (37 lbs.) medium-back
Std. seat dimensions	46 cm x 46 cm (18" x 18")
Warranty	Lifetime warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive train 18-month warranty on PG Drive control system

COLOUR CHOICES



Red



Blue

FEATURES

- 40 amp, PG GC 2 controller
- Built with ease of service in mind
- Easy side access to batteries
- Tool-less front panel access to electronics to minimize service time
- Active-Trac® Suspension

OPTIONS

- Depth-adjustable seat provides fast, simple adjustment for enhanced, personalised comfort*
- High-back seat with headrest
- High-back seat with solid seat pan
- Cup holder
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Power seat lift option

Jazzy® Select
shown in Blue



* Special order seats or non-standard dimensions will not be fold-flat and depth adjustable.

Pride
Mobility Products Ltd.

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL
(Tel) 01869 324 600 • (Fax) 01869 323 070 • www.pride-mobility.co.uk



The information contained herein is correct at the time of publication; due to the manufacturers commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.