

Smile, it's a



Parameter	Little Gem ²
Length	98cm (38.5")
Rear width	51cm (20")
Max. user weight	113 kg (250 lb)
Battery capacity	12 AH
Max safe slope	9° at 113kg (250lbs)
Turn radius	115cm (45.3")
Seat - base height	40.5 - 49.5 cm (16" - 19.5")
Max speed	6.4 kph (4 mph)
Wheel diameter front	17.5 cm (7")
Wheel diameter rear	19.0 cm (7.5")
Tyres	Solid
Charger off-board	1.5 Amps
Range	16 km (10 mls)
Overall weight	40 kg (88.1 lbs)
Seat weight	10.4 kg (22.9 lbs)
Battery weight (each)	4.31 kg (9.5 lbs)
Drive unit weight	10.7 kg (23.6 lbs)
Heaviest part (front section)	16.8 kg (37 lbs)
Ground clearance	6 cm (2.5")
Controller (pgdt)	S - drive (programmable by dealer)



Let the Little Gem²
make you smile!



Let the Little Gem²
make you smile!



Sunrise Medical Limited
 High Street, Wollaston, Stourbridge
 West Midlands DY8 4PS
 England
 Telephone: +44 (0) 1384 446622
 Fax: +44 (0) 1384 446644
 Website: www.sunrisemedical.com

INNOVATE.[™]
 BUILD FOR THE FUTURE



So easy...

To get your Little Gem² ready for transporting, simply lift off the seat, lift out the battery and lift the rear handle releasing the *Sterling Lock* that quickly separates the front and rear of the scooter. It really is that simple! (And by the way, it goes back together just as smoothly!)

The Little Gem² has been designed with your comfort and safety in mind – with *Sterling Lift* the conveniently located handles make lifting and handling the parts easier than ever before.

So comfortable...

The *Sterling Comfort* package features a soft padded seat and generous legroom; we know you'll never want to travel without your Little Gem² again!

So accessible...

With an infinitely adjustable tiller that can be pushed up and away from you, flip-up armrests and a swivelling seat, you can get on and off your Little Gem² with ease. Your Little Gem² is so compact and manoeuvrable you can take it almost anywhere.

